Detaseal® antilock

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Detaseal® antilock

GB
instructions for use

VPS precision impression silicone, highly elastic, very easy flowing

Thin flowing, very soft, elastic special silicone. Facilitates reliable demoulding and removal from the mouth in case of pronounced retention areas in the impression zone. Can be used in conjunction with VPS silicones in all undercut areas, for loose teeth with reduced stability, wide interdental spaces, or pontic bridges and bar constructions. The low hardness and high elasticity facilitate managing of even difficult areas at risk of splinting and easier demoulding.

Fig. 1

Fig. 2


Fig. 3

Fig. 4

Fig. 5

1. Preparation of the impression tray

Depending on the impression technique, select a suitable impression tray (individual or prefabricated trays made of metal or hard plastic). To avoid a push through of the impression, it may prove useful to fit stops to the tray in the non-relevant areas of the impression with wax or light-curing resin. To provide optimal bonding, we recommend coating all trays with a thin layer of **trayloc® A** (Fig. 1). Allow to dry for approx. 3 min.

2. Mixing and dosing

The Automix2-System gun is used for squeezing material from the cartridge (Fig. 2). After placing the cartridge in the mixing gun, unscrew the cartridge cap to remove. Squeeze out a small amount of material until silicone is transported evenly from both outlets (Fig. 3). To position the mixing cannula, use the guides on the cannula and the cartridge. Turn in opposite direction to lock (Fig. 4). Dosing can now be performed individually. After use, leave the mixing cannula on the cartridge until it is next used.

3. Procedure

Detaseal® antilock can be used both for correction and mixing techniques. Clean the relevant areas and dry as usual.

Double mixing technique: apply **Detaseal® antilock** to the retentive areas prior to impression taking (Fig. 5) and allow to set. Check for possible faults before proceeding to the next step. Clean, dry and then perform double mixing impression as usual.

Detaseal® antilock can also be applied together with not yet set heavy und monophasic material during the same step.

Correction technique: Here, **Detaseal® antilock** can be used in the pre-impression. To this purpose, apply the material to the undercut areas (Fig. 5) and allow to set. Prior to correction impression, check whether trouble-free repositioning of the impression is possible. If necessary, remove these parts. Clean and dry and then perform the pre-impression with putty. Avoid material overlap in the subsequent correction impression.

4. Disinfection

Rinse the impression under running, lukewarm water after removal from the mouth. Subsequent disinfection can be performed for 15 min. in 2 % glutaraldehyde or suitable, commercially available disinfectant solutions.

5. Model casting

The impression should not be cast before 30 min. Dimensional stability is guaranteed for 14 days. The recommended materials for modelling are dental gypsum of Classes III and IV as well as commercially available modelling plastics.

6. Cleaning of the tray

Remove set material mechanically with a blunt instrument. **trayloc® A** dissolves when placed in acetone or isopropanol. Only use in well ventilated rooms. Then clean and disinfect tray as usual.

Important notes on processing

- Does not combine with polyether impression compounds, condensation-linked silicones and alginates.
- Observe adequate material thickness (3 – 5 mm) when applying around loosened teeth.
- Do not use in the area of preparations.
- Solutions employed prior to impression taking (i.e. retraction solutions) may affect the setting reaction and therefore need to be removed completely.
- Latex gloves and latex-contaminated surfaces can affect the curing of **Detaseal® antilock**.
- Do not leave remains of the impression material in the patient's mouth.
- Linked impression compounds are chemically stable – avoid stains on clothing.

Miscellaneous information:

Silicone impression materials have proven themselves millions of times, adverse effects are not to be expected if used as instructed. Immune reactions such as allergies or irritations can, however, not be excluded in principle. If in doubt, we recommend to perform an allergy test prior to use.

To be used only by dental professionals.

60 sec.	1 min. 20 sec.
Processing time	Dwell time in mouth
2 min. 20 sec.	
Setting time	

Areas of application:

In undercut areas of fixed prosthetic and orthodontic restorations and around periodontally damaged loose teeth and greatly expanded interdental areas.

Technical specifications:

- Mixing volume:** 50 ml cartridges
- Dosage:** 1:1
- Product colour:** Base: pink
Catalyst: white
- Mixing time:** not applicable (Automix2-System)
- Processing time:** 60 sec.*
- Setting time:** approx. 2 min. 20 sec.*
- Dwell time in mouth:** approx. 1 min. 20 sec.
- Deformation under pressure:** approx. 22.5 %
- Recovery after deformation:** > 99.8 %
- Linear dimensional change:** approx. < 0.1 %
- Shore A hardness:** approx. 10
- Processing:** at 23°C ± 2°C, 50 ± 5 % rel. humidity
- Storage:**

* as of beginning of mixing T 23°C ± 2°C, 50 ± 5 % rel. humidity. Higher temperatures shorten the times given, colder temperatures prolong them.

Ordering information:

Detaseal® antilock
Standard pack **03535**
1 x 50 ml cartridge
6 Mixing cannulas, yellow
6 Intra-oral tips, yellow

Detaseal® hydroflow heavy
Standard pack 5:1 **02729**
Maxi-cartridge à 300 ml
Base + 60 ml catalyst
10 Dynamic mixers, yellow
1 Rotating ring, yellow

Detaseal® hydroflow putty
Jumbo pack 5:1 **03283**
Maxi-cartridge à 300 ml
Base + 60 ml catalyst
10 Dynamic mixers, yellow
1 Rotating ring, yellow

Standard pack 1:1 **02727**
2 x 250 ml jars
2 Dosing scoops

Detaseal® hydroflow mono
Jumbo pack 5:1 **03373**
Maxi-cartridge à 300 ml
Base + 60 ml catalyst
10 Dynamic mixers, yellow
1 Rotating ring, yellow

Standard pack **02754**
2 x 50 ml cartridges
6 Mixing cannulas, green

Multipack **02757**
4 x 50 ml cartridges
12 Mixing cannulas, green

Dynamic mixers yellow, 40 pieces **02833**

Rotating ring yellow, 1 piece **02836**

Mixing cannulas green, 48 pcs. **02705**

monoprint supra Standard pack **02525**
2 x 50 ml cartridges
6 Mixing cannulas, green

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